



Measuring In-vivo and In-situ Ex-vivo the 3D Deformation of the Lamina Cribrosa Microstructure under Elevated Intraocular Pressure

Junchao Wei ^a, Bin Yang ^a, Andrew P. Voorhees ^a, Huong Tran ^{a,b}, Bryn Brazile ^a, Bo Wang ^{a,b},
Joel Schuman ^c, Matthew A. Smith ^{a,b}, Gadi Wollstein ^c, and Ian A. Siga ^{a,b}




^aDepartment of Ophthalmology, University of Pittsburgh, Pittsburgh, PA, USA;

^bDepartment of Bioengineering, University of Pittsburgh, Pittsburgh, PA, USA;


^cNYU Langone Eye Center, NYU School of Medicine, New York, NY, USA.



Introduction




Introduction
Motivation
Objective
Challenges




Experimental Setup

Deformation Validations




Motivation

- Previously, ex-vivo LC deformation studies
- OCT, In-vivo imaging
- Differences between measurements



Introduction
Motivation
Objective
Challenges

Digital Volume Correlation




Objective


- To present a new in-vivo method
- To measure deformations
- To compare in-vivo and in-situ ex-vivo deformations

Introduction
Motivation
Objective
Challenges


OCT 3D Image Processing



LC Shear Strain



LC Effective Strain




$$\epsilon_{effective} = \sqrt{\frac{2}{3} \epsilon_{ij} \epsilon_{ij}}$$

$$\epsilon_{ij}' = \epsilon_{ij} - \bar{\epsilon}_{ij}': Deviatoric Strain$$

$$\bar{\epsilon}_{ij} = \frac{1}{3} \delta_{ij} \epsilon_{kk}: Hydrostatic Strain$$


Challenges in OCT

- Noise
- Lamina visibility reduced by blood vessel "shadows"
- Eye movements introduce rigid-body motions



Introduction
Motivation
Objective
Challenges

3D Image Alignment



Summary and Conclusions

- We developed a novel 3D OCT experimental and analysis technique to measure 3D LC deformations.
- The LC deforms when IOP is raised, and in-vivo deformation of the LC was larger than the in-situ ex-vivo deformation between 10 and 15mmHg.
- For studying the LC deformation contribute to glaucoma and vision loss, the ability to measure the in-vivo deformation of the LC is important.
- In-vivo measurements of LC deformation could be used as a tool for the early detection, diagnosis, and tracking of glaucoma.